

ACTIVE - TECHNICAL SCHEME ACTIVE - SPECIFICHE TECNICHE

ACTIVE - AT1A





mountable on widths from 340mm to 380mm montabile su larghezze da 340 mm a 380 mm

ACTIVE - AT2A





mountable on widths from 340mm to 380mm montabile su larghezze da 340 mm a 380 mm

ACTIVE - AT3A





mountable on widths from 340mm to 380mm montabile su larghezze da 340 mm a 380 mm

ACTIVE - AT4A





mountable on widths from 400mm to 460mm

ACTIVE - AT5A





mountable on widths from 400mm to 460mm

ACTIVE - AT6A





mountable on widths from 400mm to 460mm

Weight ± 1000gr

Supplied with 30 mm high-density foam padding and upholstery in fireproof and breathable fabric. The backrest can be mounted to the wheelchair either centrally or laterally. Both clamps guarantee a wide range of adjustments.

Peso ± 1000gr

Viene fornita un'imbottitura con schiuma ad alta densità di spessore 30mm e il relativo rivestimento in materiale ignifugo e traspirante.

Lo schienale può essere infatti montato alla carrozzina centralmente oppure lateralmente. Entrambi gli attacchi garantiscono ampie possibilità di regolazion

37.00

CLACE



adaptable to Ø 19/22/25mm tubes adattabile su tubi Ø 19/22/25mm

CLALA




adaptable to Ø 22/25mm tubes adattabile su tubi Ø 22/25mm

ARIA WHEELS SRL



DEEP - TECHNICAL SCHEME DEEP - SPECIFICHE TECNICHE


DEEP - AT1D





mountable on widths from 340mm to 380mm montabile su larghezze da 340 mm a 380 mm

DEEP - AT2D





mountable on widths from 340mm to 380mm montabile su larghezze da 340 mm a 380 mm

DEEP - AT3D





mountable on widths from 340mm to 380mm montabile su larghezze da 340 mm a 380 mm

DEEP - AT4D





mountable on widths from 400mm to 460mm montabile su larghezze da 400 mm a 460 mm

DEEP - AT5D





mountable on widths from 400mm to 460mm

DEEP - AT6D





mountable on widths from 400mm to 460mm montabile su larghezze da 400 mm a 460 mm

Weight ± 700gr

Supplied with 30 mm high-density foam padding and upholstery in fireproof and breathable fabric. The backrest can be mounted to the wheelchair either centrally or laterally. Both clamps guarantee a wide range of adjustments.

Peso ± 700gi

Viene fornita un'imbottitura con schiuma ad alta densità di spessore 30mm e il relativo rivestimento in materiale ignifugo e traspirante.

Lo schienale può essere infatti montato alla carrozzina centralmente oppure lateralmente. Entrambi gli attacchi garantiscono ampie possibilità di regolazion.

37 mm

CLACE



adaptable to Ø 19/22/25mm tubes adattabile su tubi Ø 19/22/25mm

CLALA



adaptable to Ø 22/25mm tubes adattabile su tubi Ø 22/25mm

ARIA WHEELS SRL