

BKRTA

ACTIVE - TECHNICAL SCHEME

ACTIVE - AT1A



mountable on widths from 340mm to 380mm

ACTIVE - AT2A



mountable on widths from 340mm to 380mm

ACTIVE - AT3A



mountable on widths from 340mm to 380mm

ACTIVE - AT4A



mountable on widths from 400mm to 460mm

ACTIVE - AT5A



mountable on widths from 400mm to 460mm

ACTIVE - AT6A




mountable on widths from 400mm to 460mm

Weight \pm 1000gr

Supplied with 30 mm high-density foam padding and upholstery in fireproof and breathable fabric.

The backrest can be mounted to the wheelchair either centrally or laterally.
Both clamps guarantee a wide range of adjustments.

CLACE



adaptable to \varnothing 19/22/25mm tubes

CLALA



adaptable to \varnothing 22/25mm tubes

ARIA WHEELS SRL

Registered Office_Via della Divisione Torino, 92 00143 Roma - Italy

Production Dept_Via A.Volta, 7/G 36030 Costabissara (VI) - Italy


VAT IT 12848301003 | tel.+39 0444701414

www.ariawheels.com | info@ariawheels.com

BKRTA


DEEP - TECHNICAL SCHEME

DEEP - AT1D




mountable on widths from 340mm to 380mm

DEEP - AT2D




mountable on widths from 340mm to 380mm

DEEP - AT3D




mountable on widths from 340mm to 380mm

DEEP - AT4D




mountable on widths from 400mm to 460mm

DEEP - AT5D



mountable on widths from 400mm to 460mm

DEEP - AT6D




mountable on widths from 400mm to 460mm

Weight ± 1000gr

Supplied with 30 mm high-density foam padding and upholstery in fireproof and breathable fabric.

The backrest can be mounted to the wheelchair either centrally or laterally. Both clamps guarantee a wide range of adjustments.

CLACE



adaptable to Ø 19/22/25mm tubes

CLALA



adaptable to Ø 22/25mm tubes

ARIA WHEELS SRL

Registered Office_Via della Divisione Torino, 92 00143 Roma - Italy

Production Dept_Via A.Volta, 7/G 36030 Costabissara (VI) - Italy

VAT IT 12848301003 | tel.+39 0444701414

www.ariawheels.com | info@ariawheels.com